
FebriDx[®]

FebriDx[®]

EN	English	3
FR	Français	11
NL	Nederlands	19
IT	Italiano	27
DE	Deutsche	35
ES	Español	43
PT	Português	51
EL	Ελληνικά	59

INTENDED USE

FebriDx[®] is a rapid immunoassay for the visual, qualitative, in vitro detection of elevated levels of both Myxovirus resistance protein A (MxA) and C-reactive protein (CRP) directly from fingerstick whole blood. The test measures an immune response to a suspected bacterial or viral infection in patients 1 year of age or older that present within 7 days of new onset respiratory symptoms consistent with a community-acquired acute respiratory infection.

The FebriDx[®] test aids in the diagnosis and differentiation of bacterial and viral acute respiratory infection. FebriDx[®] is not intended to diagnose any specific bacteria or virus.

The test is intended for professional use and should be used in conjunction with other clinical evidence including laboratory, radiographic, and epidemiological information.

Negative results do not preclude respiratory infection and should not be used as the sole basis for diagnosis, treatment, or other clinical and patient management decisions. In addition to utilizing radiography and clinical presentation to aid in diagnosis, additional laboratory testing (e.g. bacterial and viral culture, immunofluorescence, and polymerase chain reaction [PCR]) may be used to confirm the presence of a specific respiratory pathogen.

SUMMARY & EXPLANATION:

Acute respiratory infections (ARIs) including sinusitis, pharyngitis, bronchitis, influenza affect 20% of the population annually. The significant overlap in symptoms and signs makes it challenging for physicians to differentiate viral from bacterial infection and to identify which patients require antibiotic therapy. The vast majority of ARIs are caused by viruses, for which antibiotics provide no clinical benefit, however 30-80% receive antibiotics.¹ The over prescription of antibiotics for ARI is a leading contributor to the global antimicrobial resistance (AMR) crisis which currently causes 700K deaths annually.² FebriDx[®] utilizes dual biomarker technology to deliver high sensitivity and specificity to differentiate a viral from bacterial ARI.

The ability to rapidly identify viral infections to enable effective cohorting and prevent the spread of infection is a critical component of infection prevention control. Through differentiation of viral positive and viral negative patients, FebriDx[®] is able to optimise patient isolation decisions and enable accurate cohorting.

BIOMARKERS

MxA (Myxovirus resistance protein A)

MxA becomes elevated in the presence of acute viral infection. MxA has a low basal concentration of less than 15 ng/mL, a fast induction time of 1-2 hours, and a long half-life of 2.3 days, making it an ideal marker for viral infection.⁵ Numerous clinical studies demonstrate that MxA protein expression in peripheral blood has been shown to be a sensitive and specific marker for acute viral infection.⁵⁻¹⁰

CRP (C-reactive protein)

CRP is a nonspecific, acute-phase protein that increases during an inflammatory process, especially following severe infection. Bacterial infection is a potent stimulus of marked CRP elevation, which occurs within 4-6 hours of infection and peaks after 36 hours.^{11,12} Some viral infections, including Influenza, Adenovirus and SARS-CoV-2 cause CRP to significantly elevate.^{3,13-15}

Multiplexed Pattern of Results

In isolation, neither MxA nor CRP alone is sensitive or specific enough to differentiate viral from bacterial infection. At low levels, CRP is very sensitive but non-specific at confirming a bacterial infection. At high levels, CRP becomes very specific for bacterial infection but has low sensitivity. MxA is specific for viral infection only and is not elevated in a bacterial infection. The FebriDx[®] test produces a multiplexed pattern of results by simultaneously detecting elevated levels of MxA and CRP together to differentiate viral from bacterial infection.^{4,16-18}

PRINCIPLES OF THE TEST

The FebriDx® test is a 10-minute lateral flow immunoassay within a plastic housing that incorporates a built-in retractable lancet, blood collection and transfer tube, and buffer release mechanism. FebriDx® utilizes monoclonal anti-MxA and anti-CRP antibodies to simultaneously detect MxA at the medical decision point of approximately 40 ng/mL and CRP of approximately 20 mg/L serum equivalent.

If the fingerstick blood samples contain elevated levels of MxA or CRP at or above their respective cut-off levels, the appropriate test line will appear in the Result Window. FebriDx® also contains a control line to indicate correct sample flow and valid results.

Materials Provided

- 25 single use tests
- 1 package insert

Materials Not Provided

- Timer
- Alcohol
- Gauze
- Sterile dressing

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Keep the FebriDx® test in the sealed foil pouch until just before use. If the foil pouch is damaged do not use the test.
3. Do not use the FebriDx® test past the expiration date.
4. Use standard precautions for collecting and handling a fingerstick blood sample.
5. All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
6. Wash hands before and after performing the test and wear disposable gloves while handling specimens.
7. The lancet is sterile until the protective tab is removed. Do not use the lancet if the protective tab is not secured in place.
8. The FebriDx® test is designed to proceed in sequential order and locking mechanisms exist to prevent skipping the prior step.
9. The FebriDx® test is a single-use item with no reusable components. Proper handling and disposal methods should be established according to local, state, and federal regulations.
10. The FebriDx® test requires a visual readout. Do not interpret the test result if you have color-impaired vision.
11. A brightly lit environment is recommended for interpreting the test results.


STORAGE AND STABILITY

Store the FebriDx® test between 4-25°C (39-77°F). Unopened, the FebriDx® tests are stable until the expiration dates printed on their packaging.

ALTITUDE AND RELATIVE HUMIDITY (RH)

The FebriDx® test performed acceptably when tested at altitudes between 0-2000 meters and 5-85% RH.


TEST COMPONENTS

Unused Test

An **unused test**, or test that has not yet been activated by pressing the Buffer Release Button, will show three faint orange lines in the Result Window.

TEST PROCEDURE - Check the expiration date on all packaging.

1 Lance Finger

2 Collect Blood

3 Deliver Blood

4 Deliver Buffer

5 Check Blood Transfer

1. Tear open the foil pouch at the indicated perforation and remove the test just prior to testing.
2. Use standard precautions for collecting and handling a fingerstick blood sample.
3. Twist and pull to remove the Protective Lancet Tab. Firmly press the lancet to puncture the skin. Wipe away the first drop of blood with gauze. Gently massage finger to encourage blood flow to obtain a drop of blood that hangs from the finger.
4. Place Blood Collection Tube at a 45 degree angle below the finger, making contact with the hanging drop of blood. Avoid direct contact with the patient's finger. **ENSURE BLOOD COLLECTION TUBE IS COMPLETELY FULL.**

Note: Capillary action will automatically draw the blood sample into the Blood Collection Tube.

5. Once the Blood Collection Tube is completely filled with blood, rotate it over the Blood Transfer Zone. **ENSURE THAT BLOOD IS TRANSFERRED TO THE TEST STRIP.** Wait for ~5-10 seconds before proceeding to step #6.

Note: If the blood does not immediately transfer onto the test strip, reverse the Blood Collection Tube's rotation back to its original position and add additional blood to fill the tube. There will be some resistance. Once the Blood Collection Tube is completely filled, rotate it over the Blood Transfer Zone. **ENSURE THAT BLOOD IS TRANSFERRED TO THE TEST STRIP.** Wait for ~5-10 seconds before proceeding to step #6.

6. Lay the test on a flat surface. **Blood must be visible on the test strip before releasing buffer.** Activate the test by firmly and fully pressing the Buffer Release Button to deliver the buffer. The Buffer Release Button should be pressed within 1 minute of transferring the blood sample.

Note: If no fluid is visible within 25-30 seconds, firmly re-press the Buffer Release Button.


7. Read results at the 10-minute mark. If the blood has not cleared at 10 minutes, **wait for the background to clear before reading the test results.** Do not read the test results after 1 hour. Dispose of test in the proper biohazard receptacle.

BLOOD TRANSFER CHECK

Confirm the **Blood Transfer Confirmation Window** changed from **white** to **pink/red** to indicate the correct transfer of the blood sample.

If the **Blood Transfer Confirmation Window** does not change to a **pink/red** color, then the sample was not transferred correctly and the test needs repeating.

Blood Transfer Confirmation Window

TEST RESULTS

A **blue** control line must appear in the **Result Window** for the test to be valid.

Positive Result

The positive result lines should appear as **red** or **black** lines in the **Result Window**. An uneven or incomplete result line is due to an uneven sample distribution on the test strip. Even if the result line is faint in color, incomplete over the width of the test strip, or uneven in color, it should be interpreted as positive. A positive result indicates the presence of elevated MxA and/or CRP proteins.

* Cannot preclude co-infection. The incidence of co-infection is low and varies on severity of illness and setting.

Negative Result

If only a **blue** control line is visible in the **Result Window**, the test is deemed negative. A negative result indicates a lack of elevated MxA and CRP proteins.

Invalid Result

The absence of the **blue** control line indicates an invalid result. If an invalid result occurs, the test must be discarded and the patient retested using a new FebrIDx® test. Choose an alternative puncture site on a different finger when retesting the patient.

If the background of the **Result Window** has not cleared sufficiently for interpretation of results after 30 minutes, discard the test and retest the patient with a new FebrIDx® test.

QUALITY CONTROL

Control Line

A **blue** control line must appear in the **Result Window** for the test to be valid. The absence of a **blue** control line indicates an invalid result.

Blood Transfer Confirmation Window

Confirm the **Blood Transfer Confirmation Window** changed from white to **pink/red** to indicate the correct transfer of the blood sample.

If the **Blood Transfer Confirmation Window** does not change to a **pink/red** color, then the sample was not transferred correctly and the test needs repeating.

External Controls

External controls may be used to demonstrate that the reagents and assay perform properly. FebrIDx® external controls are available directly through Lumos Diagnostics and consist of one (1) positive control and one (1) negative control. Refer to the FebrIDx® external controls package insert for instructions on how to use. If the FebrIDx® external controls do not perform as expected, please repeat the test. If external controls fail on repeat testing do not perform patient testing and contact Lumos Diagnostics Technical Support.

LIMITATIONS

1. The FebrIDx® test is best used within three (3) days from onset of a new fever and seven (7) days from onset of new respiratory symptoms. Symptomatic patients with suspected COVID-19 may be tested up to 14 days from the onset of symptoms.
2. Only fresh capillary blood (fingerstick) must be used on the FebrIDx® test. Venous blood CANNOT be used.
3. The blood collection tube must be filled completely and applied to the test strip in order for the test to run properly. An erroneous result may occur if an insufficient blood sample is applied to the test.
4. The following conditions may lead to erroneous results:
 - Current immunosuppressive state or use of immunosuppressive drugs
 - Current use of oral anti-infective drugs
 - Current use of interferon therapy (e.g. for multiple sclerosis, HIV, hepatitis B/C)
 - Live viral immunization within the last 30 days
 - Major trauma, major surgical intervention, and severe burns within the preceding 30 days
 - Chronic fevers lasting more than 7 days that are not suspected to be associated with SARS-CoV-2
5. FebrIDx® will not identify bacterial colonization, localized infections, or periodic viral shedding without an associated systemic host response.
6. Rheumatoid Factor (RF) ≥ 100 IU/mL (normal RF: 15 IU/mL) can produce a MxA line in very rare cases.
7. Reading results before 10 minutes or after 1 hour may produce erroneous results.

EXPECTED VALUES

The prevalence of ARI varies during the year and from region to region, with outbreaks typically occurring during fall and winter. ARIs are the leading cause of morbidity, accounting for 20% of medical consultations, 30% of absenteeism, and 75% of all antibiotic prescriptions.¹⁷

Nearly half of patients presenting to the outpatient setting with respiratory symptoms do not have infection.^{16,17} A negative FebrIDx® test in a symptomatic patient can be indicative of a non-infectious illness e.g., environmental trigger of respiratory illness, allergic rhinitis, autoimmune, subclinical, or past infection. Symptomatic patients who had a negative FebrIDx® result did not have evidence of infection with standard laboratory testing^{16,17} nor returned to a healthcare provider with worsening illness that required treatment.¹⁹

CLINICAL EVIDENCE

FebrIDx® has been evaluated in multiple prospective, multicenter, blinded clinical trials with untrained operators to determine the diagnostic performance characteristics of the FebrIDx® test to identify a host immune response and differentiate viral or bacterial community-acquired febrile ARI as compared to the reference standard (standardized microbiologic and laboratory testing adjudicated by clinical experts).^{10,16,17} Subjects 1 year of age and older who presented to primary care, urgent care, or an emergency department (ED) within 3 days of an acute onset fever and within 7 days of new onset respiratory symptoms consistent with a community-acquired ARI were eligible for inclusion.

Viral testing:

- FilmArray® PCR: Influenza A/B, Adenovirus, RSV, Parainfluenza virus 1-4, Metapneumovirus, non-SARS-CoV-2 Coronavirus, and Rhinovirus
- Supplemental real-time reverse transcriptase PCR for EBV, HSV, and CMV
- EBV IgM Serology

Bacterial testing:

- FilmArray® PCR for atypical bacteria: *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae*, *Bordetella pertussis*, *Fusobacterium necrophorum*, *Neisseria gonorrhoeaea*.
- Oropharyngeal cultures (blood, chocolate, and MacConkey plates)

Laboratory testing:

- Procalcitonin (PCT) and white blood cell count (WBC), lymphocytes and percentage of immature WBC (bands)

A composite reference-testing algorithm adjudicated by an expert physician panel served as the reference standard from which FebriDx® was compared. Each patient underwent the following reference tests: (1) throat swab bacterial culture; (2) combined nasopharyngeal and oropharyngeal (NP/OP) swabs for multiplex PCR using the FilmArray® Respiratory Panel (Biomerieux, Inc.; Marcy-L'Etoile, France); (3) combined NP/OP swabs for real-time reverse transcriptase PCR for EBV, HSV, and CMV; (4) EBV IgM serum antibody with the Immunosimplicity® IS-EBV-VCA IgM test kit (Diamedix Co; Miami Lakes, FL, USA); (5) serum PCT concentration measurement using the BRAHMS PCT Kryptor™ (ThermoFisher Scientific; Waltham, MA, USA); (6) WBC with band differential, and (7) MxA protein ELISA and CRP enzyme immunoassay (Biocheck; Foster City, CA, USA). Reference testing was completed at a central laboratory and blinded to patients, treating clinicians, and study personnel who performed FebriDx® testing.

The reference testing algorithm classified patients as having a bacterial infection if any of the following 5 criteria were met: (1) throat culture positive for a bacteria that commonly causes pharyngitis (group A and C B-hemolytic Streptococci, *N. gonorrhoeaea*, *C. diphtheriae*, *A. haemolyticum*) plus PCT ≥0.1 ng/mL; (2) throat culture positive for any other bacteria plus PCT ≥0.15 ng/mL; (3) NP/OP sample PCR positive for atypical bacteria (*M. pneumoniae*, *C. pneumoniae*, *B. pertussis*) plus PCT ≥0.1 ng/mL; (4) PCT ≥0.25 ng/mL plus no identified pathogen; (5) PCT ≥0.15 ng/mL plus WBC ≥15,000 cells/mcl or presence of WBC bands plus no identified pathogen.

Pharyngeal bacterial colonization was differentiated from true systemic bacterial infection if cell culture growth occurred in the absence of an elevated PCT level (measure of host immune response). Patients with a negative FebriDx® result without an identified pathogen and a normal PCT (absent host immune response) were considered negative for infection.

The reference testing algorithm classified patients as having a viral infection if any of the following 3 criteria were met: (1) NP/OP sample PCR positive for Influenza A or B, Adenovirus, RSV, Human Metapneumovirus, Parainfluenza viruses 1-4, CMV, and HSV; (2) NP/OP sample PCR positive for EBV plus serum IgM positive for EBV; (3) PCT between 0.15 ng/mL and 0.25 ng/mL plus WBC <15,000 cells/mcl plus no WBC bands plus no identified pathogen.

Patients who did not meet the criteria for bacterial or viral infection were classified as negative by the reference testing algorithm.

FEBRIDx® PERFORMANCE DATA

Prospective, Multi-Center Clinical Studies (FebriDx® Diagnostic Performance: Bacterial vs. Viral ARI)

• ARI Cohort^{16,17}

- 429 ARI subjects symptomatic within 7 days and febrile within 3 days of presentation and who were ≥ 1 year were enrolled in the outpatient setting. 425 symptomatic ARI were included in the analyses.
- 4 subjects were excluded prior to analysis (3 had insufficient reference standard testing to determine final diagnosis and 1 had an invalid FebriDx® test)
- Results
 - 16% (66/425) Bacterial
 - 46% (196/425) Viral
 - 38% (163/425) Negative

Summary of FebriDx® Diagnostic Performance (Bacterial vs. Viral ARI)						
Study (Sample Size)	Fever (Exhibited or Reported)	Diagnosis	PPA [95% CI]	NPA [95% CI]	PPV [95% CI]	NPV [95% CI]
Shapiro (n = 121/220) ¹⁶	Exhibited on Enrollment (55%)	Bacterial	95% [77-100]	94% [88-98]	76% [59-87]	99% [93-100]
Shapiro (n = 220) ¹⁶	Reported within 3 days	Bacterial	85% [69-95]	93% [89-96]	69% [56-79]	97% [94-99]
Self (n = 205) ¹⁷	Reported within 3 days	Bacterial	80% [59-93]	93% [90-97]	63% [45-79]	97% [94-99]

PPA, Positive Percent Agreement; NPA, Negative Percent Agreement; PPV, Positive Predictive Value; NPV, Negative Predictive Value; CI, Confidence Interval

Study (Sample Size)	Fever (Exhibited or Reported)	Diagnosis	PPA [95% CI]	NPA [95% CI]	PPV [95% CI]	NPV [95% CI]
Shapiro (n = 121/220) ¹⁶	Exhibited on Enrollment (55%)	Viral	90% [81-96]	78% [62-89]	89% [82-93]	80% [67-89]
Shapiro (n = 220) ¹⁶	Reported within 3 days	Viral	90% [83-94]	76% [66-84]	83% [77-87]	85% [77-90]
Self (n = 205) ¹⁷	Reported within 3 days	Viral	87% [75-95]	83% [77-89]	64% [53-75]	95% [90-98]

PPA, Positive Percent Agreement; NPA, Negative Percent Agreement; PPV, Positive Predictive Value; NPV, Negative Predictive Value; CI, Confidence Interval

PROSPECTIVE CLINICAL STUDIES FEBRIDX® DIAGNOSTIC PERFORMANCE:

Viral infection in patients with suspected COVID-19

Two real-world diagnostic accuracy studies evaluating FebrIDx® in hospitalized subjects with suspected COVID-19 were performed during the first wave (March-April 2020) of the COVID-19 pandemic.^{3,4} Measures of diagnostic accuracy were calculated based on FebrIDx® results compared to the reference standard of either (i) SARS-CoV-2 PCR on combined nose and throat swabs³ and/or SARS-CoV-2 serum antibody test⁴ as well as the European Centre for Disease Prevention and Control (ECDC) and U.S. Centers for Disease Control and Prevention (CDC) probable case definition for COVID-19 (any person meeting the clinical criteria with an epidemiological link or any person meeting the diagnostic criteria [radiological evidence consistent with COVID-19]).^{20,21}

- 248 symptomatic subjects with suspected COVID-19 were enrolled in the ED and in-patient setting³
 - o 48% (118/248) of the cohort had SARS-CoV-2 confirmed by PCR (ECDC/CDC definition for confirmed COVID-19)
 - o 50% (125/248) of the cohort had met the ECDC/CDC case definition for probable (n = 7) or confirmed (n = 118) COVID-19
- 47 symptomatic subjects with suspected COVID-19 were enrolled in the ED and in-patient setting⁴
 - o 66% (31/47) had SARS-CoV-2 confirmed by PCR
 - o 74% (35/47) had a diagnosis of COVID-19 (n = 31 confirmed by PCR; n = 3 met ECDC/CDC case definition for probable COVID-19 together with antibodies)

RESULTS

	n	Sensitivity [95% CI]	Specificity [95% CI]	PPV [95% CI]	NPV [95% CI]
FebrIDx® v. SARS-CoV-2 PCR ¹¹	248	93% 110/118 [87-97]	86% 112/130 [79-92]	86% 110/128 [79-91]	93% 112/120 [87-97]
FebrIDx® v. Probable Case Definition (ECDC/CDC) ¹¹	248	93.6% 117/125 [87.8-97.2]	91.1% 112/123 [84.5-95.5]	91.4% 112/120 [84.5-95.5]	93.3% 117/128 [88.7-96.8]
SARS-CoV-2 PCR, Positive Antibody Test, or Probable Case Definition (ECDC/CDC) ¹¹	47	100% 34/34 [87.4-100]	100% 13/13 [71.7-100]	100% 34/34 [84.7-100]	100% 13/13 [71.7-100]
FebrIDx® v. Bacterial Infection ¹⁶	47	100% 8/8 [59.8-100]	92.3% 36/39 [78-98]	72.7% 8/11 [39.3-92.3]	100% 36/36 [88-100]

PRECISION AND REPRODUCIBILITY STUDY

Samples were prepared in fresh EDTA whole blood with recombinant MxA and CRP proteins. Six (6) samples, consisting of a combination of no analyte, C5 and C95 concentrations of MxA and CRP were tested.²² A total of 1080 determinations were performed by untrained operators at three (3) different sites over five (5) days. The study demonstrates overall reproducibility among three (3) lots of material, among three (3) separate sites, and among six (6) separate users.

Study	Results
Within Run Reproducibility	83.3*-100%
Within Day Reproducibility	91.7-100%
Within Laboratory Reproducibility	96.7-99.2%
Overall Reproducibility	97.8%

*One day, one operator. Otherwise within run reproducibility ranged 91.7-100%.

INTERFERING SUBSTANCES

The Interfering Substances Verification Study assessed the impact of substances that might be found in samples on the analytical specificity and analytical sensitivity of the FebrIDx® test. This assessment was performed by evaluating three replicates each of a series of samples that included MxA and CRP at the C95 concentration (i.e. low positive) and negative levels in whole blood, spiked with interfering substances. Positive and negative interference with the potentially interfering substances was evaluated by three independent researchers blinded to the sample composition.

The following substances were evaluated on the FebrIDx® test and found to not interfere at the listed test concentrations:

Test Substance	Concentration
Acetaminophen	15.6 mg/dL
Acetylsalicylic acid	3 mg/dL
Alcohol	789 mg/dL
Azithromycin	1.11 mg/dL
Biotin	3500 ng/mL
Caffeine	10 mg/dL
Celecoxib	0.879 mg/dL
Cetirizine HCl	0.435 mg/dL
Conjugated Bilirubin	40 mg/dL
Dextromethorphan	1.56 ug/dL
Doxycycline	1.8 mg/dL
Furosemide	1.59 mg/dL
HAMA	524.6 ng/mL
Hemoglobin	1000 mg/dL
Ibuprofen	21.9 mg/dL

Test Substance	Concentration
Imipenem	18 mg/dL
Levofloxacin	3.6 mg/dL
Loratadine	0.5 mg/dL
Nicotine	0.097 mg/dL
Oxymetazoline HCl	0.09 mg/dL
Phenylephrine	0.003 mg/dL
Prednisolone	0.120 mg/dL
Protein (total)	9 g/dL
Rheumatoid Factor (RF)	50 IU/mL
Salmeterol	6.03 ug/dL
Tiotropium	4.80 ng/dL
Triglycerides	1500 mg/dL
Unconjugated Bilirubin	40 mg/dL
Vancomycin	12 mg/dL

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1. Harris AM, Hicks LA, Qaseem A. Appropriate antibiotic use for acute respiratory tract infection in adults: advice for high-value care from the American College of Physicians and the Centers for Disease Control and Prevention. *Annals of Internal Medicine* 2016;164:425-34.
2. O'Neill J. Tackling drug-resistant infections globally: final report and recommendations: Review on antimicrobial resistance; 2016.
3. Clark TW, Brendish NJ, Poole S, et al. Diagnostic accuracy of the FebriDx host response point-of-care test in patients hospitalised with suspected COVID-19. *Journal of Infection* 2020.
4. Karim N, Ashraf MZ, Naeem M, et al. Utility of the FebriDx point-of-care test for rapid triage and identification of possible coronavirus disease 2019 (COVID-19). *Int J Clin Pract*. 2020.
5. Nakabayashi M, Adachi Y, Itazawa T, et al. MxA-based recognition of viral illness in febrile children by a whole blood assay. *Pediatric research* 2006;60:770-4.
6. Itazawa T, Adachi Y, Imamura H, et al. Increased lymphoid MxA expression in acute asthma exacerbation in children. *Allergy* 2001;56:895-8.
7. Forster J, Schweizer M, Schumacher RF, Kaufmehl K, Lob S. MxA protein in infants and children with respiratory tract infection. *Acta Paediatrica* (Oslo, Norway : 1992) 1996;85:163-7.
8. Engelmann I, Dubos F, Lobert P-E, et al. Diagnosis of viral infections using myxovirus resistance protein A (MxA). *Pediatrics* 2015;135:e985-e993.
9. Halminen M, Ilonen J, Julkunen I, Ruuskanen O, Simell O, Makela MJ. Expression of MxA protein in blood lymphocytes discriminates between viral and bacterial infections in febrile children. *Pediatric research* 1997;41:647-50.
10. Sambursky R, Shapiro N. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. *European clinical respiratory journal* 2015;2:28245.
11. Andreola B, Bressan S, Callegaro S, Liverani A, Plebani M, Da Dalt L. Procalcitonin and C-reactive protein as diagnostic markers of severe bacterial infections in febrile infants and children in the emergency department. *The Pediatric infectious disease journal* 2007;26:672-7.
12. Simon L, Gauvin F, Amre DK, Saint-Louis P, Lacroix J. Serum procalcitonin and C-reactive protein levels as markers of bacterial infection: a systematic review and meta-analysis. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2004;39:206-17.
13. Halminen M, Ilonen J, Julkunen I, Ruuskanen O, Simell O, Mäkelä MJ. Expression of MxA protein in blood lymphocytes discriminates between viral and bacterial infections in febrile children. *Pediatric research* 1997;41:647.
14. Salonen EM, Vaheri A. C-reactive protein in acute viral infections. *Journal of medical virology* 1981;8:161-7.
15. Chew, Keng Sheng. "What's new in Emergencies Trauma and Shock? C-reactive protein as a potential clinical biomarker for influenza infection: More questions than answers." *Journal of emergencies, trauma, and shock* 5.2 (2012):115.
16. Shapiro NI, Self WH, Rosen J, et al. A prospective, multi-centre US clinical trial to determine accuracy of FebriDx point-of-care testing for acute upper respiratory infections with and without a confirmed fever. *Annals of medicine* 2018;50:420-9.
17. Self WH, Rosen J, Sharp SC, et al. Diagnostic accuracy of FebriDx: a rapid test to detect immune responses to viral and bacterial upper respiratory infections. *Journal of clinical medicine* 2017;6:94.
18. Davidson M. FebriDx Point-of-Care Testing to Guide Antibiotic Therapy for Acute Respiratory Tract Infection in UK Primary Care: A Retrospective Outcome Analysis. *J Infect Dis Preve Med* 2017;5:165.
19. Lumos Diagnostics internal data, 2019.
20. The Centre for Disease Prevention and Control (ECDC). Case definition for coronavirus disease 2019 (COVID-19), as of 29 May 2020. Available at: <https://www.ecdc.europa.eu/en/covid-19/surveillance/case-definition2020>.
21. The Centers for Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19): 2020 Interim Case Definition, Approved April 5, 2020. Available at: <https://www.cdc.gov/nndss/conditions/coronavirus-disease-2019-covid-19/case-definition/2020/>.
22. Food and Drug Administration (FDA). Assay Migration Studies for In Vitro Diagnostic Devices: Guidance for Industry and FDA Staff. April 25, 2013:1-47. Available at: <https://www.fda.gov/medial/73669/download>.

ORDERING AND CONTACT INFORMATION / INFORMACIÓN DE CONTACTO Y PARA PEDIDOS / BESTELINFORMATIE EN CONTACTGEGEVENS / INFORMATIONS DE COMMANDE ET COORDONNÉES / BESTELL- UND KONTAKTINFORMATIONEN / INFORMAZIONI PER GLI ORDINI E CONTATTI / INFORMAÇÕES DE CONTACTO E PARA ENCOMENDAS / ΠΛΗΡΟΦΟΡΙΞ ΠΑΡΑΓΓΕΛΙΩΝ ΚΑΙ ΕΠΙΚΟΙΝΩΝΙΑΣ

Telephone: +1.855.LumosDx

Email: technical.support@lumosdiagnostics.com

Manufacturer and United States Representative / Fabricante y representante en Estados Unidos / Fabrikant en gemachtigde in de Verenigde Staten / Fabricant et mandataire pour les États-Unis / Hersteller und Bevollmächtigter in den USA / Produttore e rappresentante USA / Fabricante e representante nos EUA / Κατασκευαστής και αντιπρόσωπος για τις Ηνωμένες Πολιτείες της Αμερικής

European Representative / Representante europeo / Europees Gemachtigde / Mandataire pour l'Europe / Bevollmächtigter in der EU / Rappresentante per l'Europa / Representante na Europa / Αντιπρόσωπος για την Ευρώπη

■ Lumos Diagnostics, Inc.
2724 Loker Avenue West
Carlsbad, CA 92010 USA
t: +1.855.LumosDx
LumosDx.com

■ MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert, Germany
t: +49.6894.581020 f: +49.6894.581021
Mt-procons.com

Lumos Diagnostics is a trade name of Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. is a wholly owned subsidiary of Lumos Diagnostics Holdings LTD. / Lumos Diagnostics es un nombre comercial de Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. es una filial de titularidad exclusiva de Lumos Diagnostics Holdings LTD. / Lumos Diagnostics is een handelsnaam van Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. is een volle dochteronderneming van Lumos Diagnostics Holdings LTD. / Lumos Diagnostics est une marque commerciale de Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. est une filiale en propriété exclusive de Lumos Diagnostics Holdings LTD. / Lumos Diagnostics ist eine Handelsbezeichnung von Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. ist einem Tochterunternehmen von Lumos Diagnostics Holdings LTD. / Lumos Diagnostics è un nome commerciale di Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. è una consociata interamente controllata da Lumos Diagnostics Holdings LTD. / Lumos Diagnostics é uma marca comercial da Lumos Diagnostics Holdings LTD. Rapid Pathogen Screening, Inc. é uma subsidiária totalmente detida da Lumos Diagnostics Holdings LTD. / H Lumos Diagnostics αποτελεί εμπορική ονομασία της Lumos Diagnostics Holdings LTD. H Rapid Pathogen Screening, Inc. είναι μια εξολοκλήρου ελεγχόμενη θυγατρική της Lumos Diagnostics Holdings LTD.

©2023 Rapid Pathogen Screening, Inc. All rights reserved. / Reservados todos los derechos. / Alle rechten voorbehouden. / Tous droits réservés. / Alle Rechte vorbehalten. / Tutti i diritti riservati. / Todos os direitos reservados. / Με την επιφύλαξη παντός δικαιώματος.